



Arun Bantval <arbantval@gmail.com>

New Protocols

Jack Bouchakjian <jackb.prs@gmail.com>

Fri, Mar 20, 2020 at 2:47 AM

To: Ankit Parikh <aparikh3191@yahoo.com>

Cc: Riding <riding@plainsbororescuesquad.org>, "Reichard, Daniel" <dreichard@plainsboropolice.com>, Brian Barbarise <bbarbarise@plainsboropolice.com>, Bryan Boehle <bboehle@plainsboropolice.com>, Kenneth Mitchell <kmitchell@plainsboropolice.com>

Thanks Ankit.

Adding to the N95 masks, it is best to place the mask on your face prior to applying the straps.

You then need to pinch the metal strip to mold it against your nose/face.

To check for a good seal, you can cover the mask with both hands and breath in and out to feel vacuum/pressure.

As mentioned, removing facial hair is the right way.

Thanks,
Jack

[Quoted text hidden]